

OPHTHALMIC UNIT  
**Elegance**



## REFRACTION UNIT

**ELEGANCE** it is a modern support station for 2 instruments fully compatible with a wide range of ophthalmic devices. The refraction unit is equipped with a practical and ergonomic command console which allows the optimal use of the instrumentation and gives the possibility to memorize all the working functions of the device for a fast and comfortable use. As the whole range of CSO ophthalmic units, ELEGANCE can be completely customized with a wide choice of colors.




Eye charts projector support



Remote control compartment



Commands console




Ceiling light with three light sources suitable for every unit with:

- LED light with adjustable intensity
- with light spot
- upper ambient LED light (optional)




- Phoropter support (optional)



## STANDARD EQUIPMENT

- Central structure with command console
- Column with adjustable LED spotlight with variable intensity
- Eye charts projector support
- Plane with manual movement, 2 positions
- Retinoscope/ophthalmoscope compartment



## ACCESSORIES (OPTIONAL)

- Chair with electronic height adjustment
- Self-balancing arm for phoropter + support
- Ceiling light with 3 LED lights
- Movable LED spotlight
- Arrangement video system
- LCD monitor support
- Plane for auto refractometer
- Plane for column monitor
- Desk
- Drawers
- Version with opening on the left
- Special supports for chin rest


## ELECTRICAL DATA

|                                             | Version 230V                                                                              | version 120V       |
|---------------------------------------------|-------------------------------------------------------------------------------------------|--------------------|
| Supply voltage                              | 230 Vac $\pm 10\%$                                                                        | 120 Vac $\pm 10\%$ |
| Power frequency                             | 50Hz                                                                                      | 60Hz               |
| Network fuses                               | 4A T 5x20mm                                                                               | 8A T 5x20mm        |
| Auxiliary sockets (chair and devices)       | 230Vac 50Hz                                                                               | 120Vac 60Hz        |
| Max power absorbed by the auxiliary sockets | 100 VA                                                                                    | 100 VA             |
| Max power absorbed by the unit:             | 600 VA                                                                                    | 600 VA             |
| Voltage for the table instruments           | 6Vac or 12Vac Ophthalmometer   Slit Lamp<br>6Vac or 12Vac Fixation point of the Chin rest |                    |


# COLORS\*

They make the difference and give personality to the studio furniture.

## COLORS LACQUERED - optional

|                                                                                     |                                                                                     |                                                                                      |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
|    |    |    |
| cod. 107                                                                            | cod. 132                                                                            | cod. 135                                                                             |
|    |    |    |
| cod. 156                                                                            | cod. 162                                                                            | cod. 223                                                                             |
|    |    |    |
| cod. 218                                                                            | cod. 208                                                                            | cod. 187                                                                             |
|    |    |    |
| cod. 757/296-D                                                                      | cod. 182                                                                            | cod. 131                                                                             |
|  |  |  |
| cod. 118                                                                            | cod. 120                                                                            | cod. 130                                                                             |
|  |  |  |
| cod. 228                                                                            | cod. 236                                                                            | cod. 849                                                                             |
|  |  |  |
| cod. 945/C                                                                          | cod. 244                                                                            | cod. 259                                                                             |

## LACQUED COLORS - standard

|                                                                                     |                                                                                     |                                                                                      |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
|  |  |  |
| cod. 757/296-E                                                                      | cod. 889/F                                                                          | cod. 857/C                                                                           |

\* The colors are just approximate.  
Intensity color variations can be caused by the chemicals reactions of the used matters.

*The specifications and the images are not contractually binding and can be modified without notice.*



*YOUR PROFESSIONAL PARTNER SINCE 1967*



Via degli Stagnacci 12/E  
50018 - Scandicci - FI - Italy  
tel +39 055 72219 | fax +39 055 7215557  
email. [cso@csoitalia.it](mailto:cso@csoitalia.it) | web. [www.csoitalia.it](http://www.csoitalia.it)

